



Influence of Aggregates on *In Vitro* T Cell Responses

Christine Bryson

Antitope

EIP Tuesday 7th February

Aggregation of Protein Therapeutics

- Highly undesirable; decreased biological activity of the protein, and the potential to trigger adverse immune responses.
- Enhanced immune responses to protein aggregates have been reported in animal and clinical studies (*van Beers, Pharma. Res. 2011, Moore and Leppert J. Clin. endocrin. 1980*).
- Aggregates displaying non-native protein conformations may be seen by the immune system as neoantigens and facilitate breaking of tolerance which could trigger antibody formation.
- Characteristic adaptive humoral immune response requiring T cell activation.

Investigation into Effects of Protein Aggregates on T Cell Immune Responses

- Induce protein therapeutics to aggregate using variety of 'stress' conditions.
- Compare ability of different aggregate samples to induce in vitro T cell responses.
- Aggregates produced with different properties:
 - % aggregation in sample
 - Soluble/insoluble
 - Affect structure i.e. covalent (loss of 2° and/or 3° structure)
 - Small <1µm or large >10µm (typically hydrophobic)


Study 1:

In Vitro T Cell Activation Induced by Human IgG and Human Insulin Aggregates

- PBMC (depleted of CD8+ T cells) from 50 healthy donors stimulated with aggregate samples and T cell activation measured.
- IgG antibody induced to aggregate under stress conditions:
 - Shaking 24h at 500rpm
 - 5 cycles of freeze thawing
- Insulin aggregates induced under stress conditions:
 - Glutaraldehyde cross linking
 - Heat denaturation

Collaboration: Dr Jiskoot, University Leiden

In Vitro EpiScreen Time Course T Cell Assay: CBA, Proliferation and ELISpot




Study 1:


Overview Generation of IgG1 Aggregates

ID	Formulation	Stress conditions	Structural properties	Quantity /Type of aggregates
F/T	1 mg/mL IgG 100 mM phosphate pH 7.2	5 freeze-thawing cycles	<ul style="list-style-type: none"> • non-covalent aggr. • intact secondary & tertiary structure 	<p>< 0.1% particles larger 1μm</p> <p>0.2% larger soluble aggr.</p> <p>1.8% dimers</p> <p>97.3% monomer</p>
Control F/T		Unstressed control	native IgG	<p>1.5% dimers</p> <p>98.0% monomer</p>
Shake	1 mg/mL 10 mM citrate 5% (w/v) Sucrose pH 6	Shaking 500 rpm, 24h	<ul style="list-style-type: none"> • mainly non-covalent aggr. • loss of secondary and tertiary structure • hydrophobic aggr. 	<p>14% large insoluble aggr.</p> <p>0.5% soluble aggr.</p> <p>1.3% dimer/trimer</p> <p>84% monomer</p>
Control Shake		Unstressed control	native IgG	<p>1.5% dimers</p> <p>98.0% monomer</p>

Study 1: In Vitro T Cell Stimulation with IgG1 Aggregates



Magnitude of Response



Box Whiskers (Whiskers 10-90 percentile)

Frequency of Positive Responses (%)

SI ≥ 2.00, p < 0.05 (cpm or spw test vs baseline)

	Control F/T	F/T	Control Shake	Shake
Frequency				
Proliferation	4%	10%	6%	32%
IL-2	6%	8%	10%	24%
Both	4%	8%	6%	24%

2 fold

4 fold

Study 1: In Vitro T Cell Stimulation with Insulin Aggregates


Frequency Positive Responses

SI \geq 2.00, $p < 0.05$ (cpm or spw test vs baseline)

	Control Glu	Glu	Control Heat	Heat
Frequency				
Proliferation	4%	16%	6%	12%
IL-2	4%	16%	8%	16%
Both	2%	16%	6%	12%

8 fold

2 fold

Insulin Study  Glutaraldehyde Stress
High frequency soluble aggregates

IgG Study  Shaking Stress (IgG Study)
High frequency large insoluble aggregates

Loss of
tertiary
structure


Study 2:

Humira and Herceptin Aggregates


- Samples comprised of monomeric or aggregated clinical antibodies.
- Compare aggregated samples of antibodies associated with high and low risk of clinical immunogenicity in *in vitro* T cell assays.
- Test samples for T cell immunogenicity using EpiScreen time course T cell assay against cohort of 50 donors.
- Measure T cell activation by proliferation, IL-2 ELISpot.
- Determine relative risk of immunogenicity.

Study 2: Donor Selection


50 healthy donors




Study 2: Aggregation Induces High Magnitude T Cell Responses



Study2: T Cell Precursor Against Aggregated Antibody



Herceptin




Humira

Study 1 and 2: Conclusions

- Aggregates of antibodies and other protein therapeutics enhance in vitro T cell responses.
 - Antibodies that show different levels of clinical immunogenicity (Humira and Herceptin) have a similar capacity to activate T cells when induced to aggregate.
 - The kinetics of T cell proliferation indicate that aggregates can induce rapid T cell (polyclonal?) responses.
- Mechanism of enhanced T cell responses with aggregates, is likely to include enhanced activation of DC: Possibly via PRR, leading to increased expression of co-stimulatory molecules and proinflammatory cytokine secretion.

Characterisation of Aggregate Induced T Cell Response by Cytometric Bead Array (CBA)




1. Capture beads have a unique fluorescence intensity and are coated with a capture antibody specific for a single cytokine.
2. A combination of different beads is mixed with a sample or standard + mixture of detection antibodies that are conjugated to a reporter molecule (PE).
3. Samples are acquired on a flow cytometer.
4. Standard curves are generated and test sample concentrations determined.

Study 3: Cytokine Secretion from Aggregate Treated CD8⁺ Depleted PBMC and DC

- Cytokine secretion from aggregate and monomer stimulated PBMC and aggregate stimulated MØ-derived dendritic cells
- Characterised aggregates of a clinical antibody: Stir, freeze thaw, and heat stress.
- Included metal catalyst oxidation. Strong responses in animal models (*van Beers 2011, Hermeling 2005 Pharma. Res*) .
- Cytokine secretion was measure by multiplex cytokine bead assay.

Study 3, Cytokine Profile: Antibody Aggregate Stimulated CD8⁺ Depleted PBMC







Slide 17

CB1

Will probably take out MO. Will speak to Andrea first.

Christine Bryson, 2/3/2012


Study 3, Cytokine Profile: Antibody Aggregate Treated MØ-Derived Dendritic Cells



Cytokine Secretion Summary

- Shaking/stirred samples consistently show enhanced T cell activation – implications for transport, handling.
- Metal oxidised sample confirmed animal study results (van Beers 2011)
- Enhanced secretion of proinflammatory cytokines from DC.
- Preliminary Data suggests upregulation of CD80, CD86, CD40L.
- IL-10 upregulated in both PBMC and DC cultures. Suggesting activation of monocytes through PRR.
- DC increased IL-1B – mechanism for lymphocyte activation. Suggests a role for NLR rather than TLR?

Proposed Mechanism of Action



Conclusions

- Protein aggregates can enhance *in vitro* T cell stimulation.
- Enhanced stimulation influenced by aggregate e.g. size, quantity and/or physiochemical properties.
- Evidence of enhanced DC activation status. Enhanced uptake, potential multimeric interaction with PRR (TLR or NLR?), increased expression of co-stimulation molecules, proinflammatory cytokine secretion.
 - Hou et al 2008 – CpG multimers bind to TLR9 = enhanced DC maturations – T cell response.
- *Future directions-Tracking fluorescent labelled aggregate through APC, DC phenotyping, Caspase 1 inhibitor, TLR blocking antibodies.*

Contributors

Antitope

Matthew Baker
Tim Jones

Immunology Group

Ed Cloake
Helen Reynolds
Katie Welch

Beverley Campbell
Tania Fisher
Kirti Patel
Marina Coats
Matthew Gothard
Alyson Rust
Emma Pullen

Collaborators

University of Leiden

Wim Jiskoot
Andrea Hawe
Vasco Fillipe
Stefan Romeijn

